

## D1F05A - NOTES AND COMPUTED VARIABLES

### (Phlebotomy)

#### New Release of Version 4.1

When the samples drawn for glucose were analyzed, a few data entry errors were found in variable D05TB1E. These values have been corrected in this release (D1F05A). Also, variable labels for the individual tubes were incorrectly labeled in the previous release.

#### Extreme Fasting Values

Some values of D05FSTMN (minutes fasting) are unusually small or large. All values have been verified with the clinical centers.

#### Computed Variables

D05FSTMN: total number of minutes fasting based on time and date an individual last ate and time and date blood drawn

D05FAST: fasting status based on the calculated variable D05FSTMN. D05FAST was generated using the following statements to indicate a 12 hour fast:

```
IF 720<=D05FSTMN<=2940 THEN D05FAST = 'YES';  
IF 0 < D05FSTMN < 720 THEN D05FAST = 'NO';
```